



● Detection and semi-quantitative determination of designer benzodiazepines in serum using LC-MSⁿ

Benzodiazepines play an important role in forensic toxicology as they are widely prescribed in the treatment of psychiatric disorders and also used as drugs of abuse. In 2012 the group of New Psychoactive Substances (NPS) including numerous synthetic cannabinoids and designer stimulants (“bath salts”) was extended by the inclusion of benzodiazepine-type compounds.

Introduction

Benzodiazepines such as phenazepam and etizolam - which are still prescribed in some countries - were initially available on the Internet providing an attractive and alternative source to prescription-only benzodiazepines.

In the last years, the group of so-called designer benzodiazepines has been enlarged by compounds that are either precursors (e.g. diclazepam) or active metabolites (e.g. norfludiazepam) of known benzodiazepines or that combine structural properties of different

classical benzodiazepines (e.g. flubromazolam). Since patents and scientific literature describe the synthesis, and detail results of animal model studies, for more than a hundred different benzodiazepines, it can be assumed that this sub-group of NPS will extend quickly in the

Keywords:
Screening of benzodiazepines, amphetamine speed, semi-quantitative determination

future. In this study we describe the workflow to augment the Toxtyper™ library with new compounds of interest such as synthetic drugs. In addition, we show how semi-quantitative results may be obtained within the Toxtyper screening workflow.

Methods

Sample Preparation:

Serum samples were extracted using an alkaline liquid-liquid extraction currently used for other routine LC-MS methods in the authors' laboratory¹¹. D5-Diazepam, D4-Haloperidol and D3-Doxepin, 50 ng/mL each, were added as internal standards (IS) prior to extraction. 1 mL of serum was extracted using 0.5 mL of 1-chlorobutane and 1.5 mL of borate buffer (pH 9). After 5 min of mixing, followed by a 10 min centrifugation at RT, the organic supernatant was transferred to an LC-vial, dried under a gentle stream of nitrogen (40°C) and the residue dissolved in 25 µL eluent A : B (50 : 50 (v/v)). This sample preparation is identical to that used for the Toxtyper™ screening, thus these extracts can be re-used for either comprehensive screening or semi-quantitative analysis of benzodiazepines.

LC-MS conditions:

For screening of benzodiazepines in serum, the Toxtyper LC gradient² was modified slightly to optimize the elution over the LC runtime, and to specifically decrease co-elution of isobaric compounds or co-elution of benzodiazepines with isobaric isotopes. As all benzodiazepines ionize in positive ESI mode, zero delay polarity switching of the Toxtyper approach was turned off and the ESI source was operated in positive mode only. The adapted LC-MS conditions are shown in [Table 1](#) and [Table 2](#).

Table 1: LC-conditions for separation of benzodiazepines.

LC-Conditions	
LC-System	Dionex UltiMate 3000 LC-System
Eluent A	Water, 2 mM ammonium formate, 0.1% formic acid, 1% acetonitrile
Eluent B	Acetonitrile, 2 mM ammonium formate, 0.1% formic acid, 1% water
Analytical column	Acclaim® RSLC 120 C18 2,2 µm 120A 2.1x100 mm
Flow rate	500 µL/min
Injection volume	2 µL
Gradient:	0.0 to 0.2 min: 1% B 0.2 to 0.5 min: 1% B to 35% B, linear 0.5 to 6.0 min: 35% B to 40% B, linear 6.0 to 8.5 min: 40% B to 95% B, linear 8.5 to 11.0 min: 95% B 11.0 to 11.1 min: 95% B to 1% B, linear 11.1 to 13.0 min: 1% B

Table 2: MS-conditions for data dependent acquisition of spectra.

MS-Conditions	
MS-System	amaZon speed™ ion trap
Ion source	ESI source, positive mode
Scan mode	UltraScan: 70 - 600 Da at 32.500 Da/s
	Auto MS ⁿ mode: n = 3
	Scheduled Precursor List to trigger MS ² - and MS ³ - spectra




Figure 1: Evaluation of matrix effects in human serum (n = 5).

Spectra recording and building of the library:

For identification of new compounds of interest, spectral information from the Toxtyper library was extended by the addition of spectra of recently emerging designer benzodiazepines. Methanolic solutions of so called "research chemicals" were diluted in eluent (A : B = 50 : 50, v/v) and 1 µL was injected into the LC-MS-system

using the LC conditions described above. To acquire library spectra, the most intense m/z ratio was isolated for generation of MS² or MS³ spectra without a scheduled precursor list (SPL) and without 'active exclusion' of precursors. Within the Toxtyper workflow, the best spectrum was selected from the compound list, deconvoluted, converted to line spectrum, and checked according to the following criteria: mass difference

of the precursor ion in MS¹ and the calculated monoisotopic mass is within ± 0.25 amu; the intensity of the base peak in MS²-stage is at least 10⁵ counts, and the intensity of the base peak in MS³-stage is at least 10⁴ counts. Spectra matching these criteria were added to the library, including meta data such as empirical formula and associated structure.

Evaluation of limits of detection (LOD):

The limit of detection was evaluated by the addition of different concen-

trations of benzodiazepines to pooled blank serum aliquots (n = 6, 1 mL each) with different concentrations of benzodiazepines. The LOD of a substance was set at the lowest concentration resulting in a positive automatic identification in replicate determination (Table 3).

Evaluation of matrix effects:

Matrix effects (ME) and recovery (RE) were assessed according to Matuszewski et al.^[3] using known blank serum samples from five volunteers. For all three sets, two replicates of a

low and high concentration level were prepared and subsequently analyzed. Average ME varied between 53 and 211% (SD: 3.0 - 33.6) for the low concentration levels (25 ng/mL, 50 ng/mL for compounds with high LOD) and between 68 and 244% (SD: 1.9 - 24.4) at 250 ng/mL.

Evaluation of ME (110%, SD: 11.4) and RE (0%) for the high concentration of nitrazepam indicated that the non-satisfactory LOD is likely caused by compound degradation through contact with the serum or the extraction solvents of the LLE.

Table 3: Limits of detection (LOD) in pooled human serum.

Name	LOD [ng/mL]	Name	LOD [ng/mL]	Name	LOD [ng/mL]
2-OH-Ethylflurazepam	25	Demoxepam	50	Norflunitrazepam	5
3-OH-Bromazepam	50	Desalkylflurazepam	10	Nifoxipam	-
3-OH-Flubromazepam	25	Deschloroetizolam	1	Nimetazepam	10
3-OH-Phenazepam	25	Diazepam	5	Nitrazepam	25
7-Aminoclonazepam	10	Diclazepam	1	Nitrazolam	5
7-Aminoflunitrazepam	1	Etizolam	1	Norclobazam	10
7-Aminonitrazepam	1	Flubromazepam	5	Nordazepam	5
Adinazolam	5	Flubromazolam	5	Oxazepam	25
α-OH-Alprazolam	5	Fludiazepam	10	Phenazepam	5
α-OH-Midazolam	5	Flunitrazepam	1	Prazepam	1
α-OH-Triazolam	5	Flunitrazolam	5	Pyrazolam	5
Alprazolam	1	Flurazepam	1	RO-5-4864	10
Bromazepam	10	Fonazepam	5	Temazepam	1
Chlordiazepoxide	5	Loprazolam	1	Tetraazepam	5
Clobazam	1	Lorazepam	5	Triazolam	5
Clonazepam	5	Lormetazepam	5	Zaleplone	5
Clonazolam	5	Meclonazepam	5	Zolpidem	1
Cloniprazepam	1	Medazepam	5	Zopiclone	5
Clotiazepam	1	Metizolam	1		
Delorazepam	5	Midazolam	1		




Figure 2: MS-conditions for data dependent acquisition of spectra.

Evaluation of linear concentration range:

To establish the semi-quantitative aspect of the target screening, the linear range of each analyte must be evaluated and defined. 1 mL pooled human serum (n = 8) was spiked with 1.0 to 500 ng/mL of each compound. To assess linearity, the peak area of each molecular ion ($[M+H]^+$) was used, without normalization using the peak area of internal standards (IS), as this approach resulted in higher R^2 values. Nevertheless, the use of internal standards is crucial for semi-quantitative analysis in analytical laboratory settings.

The majority of compounds showed linear calibration curves from 5.0 ng/mL to 500 ng/mL, as exemplified above for adinazolam (Figure 3, upper panel). Desalkylflurazepam (Figure 3, central panel) is shown as an example of good linearity over the whole concentration range including one outlier, likely due to matrix interference. For 20 of the analytes studied, R^2 values > 0.99 were observed, and this number increased to 38 with the elimination of one outlier per compound.

Evaluation of semi-quantitative results in serum:

A single-point calibrator (c = 50 ng/mL) and an internal standard (D5-Diazepam) each in pooled serum, were used for semi-quantitative screening of serum samples. The peak area ratio of the molecular ion of the analyte and the internal standard was used for quantitation. Data evaluation was automatically carried out by the DataAnalysis software package. Positive findings below or above the linear range were reported as ' $< cal_{Low}$ ' or ' $> cal_{High}$ ', respectively. Following good laboratory practice, automatically generated screening results must be reviewed for infrequent false positives by manual inspection of the spectra available within the report.




Figure 3: Examples of calibration curves without normalization using an isotope labeled standard.




Figure 4: Results of semi-quant at 5 ng/mL.

Results and Discussion

The current spectral library contains 34 prescription benzodiazepines (most commonly prescribed in Germany), 21 designer benzodiazepines and three z-drugs, enabling screening for 58 benzodiazepine-type substances. The library can easily be extended with the emergence of new designer benzodiazepines entering the drug market or according to specific needs of the user. The limit of detection was 5 ng/mL for the majority of the analytes, while nine compounds could only be detected at concentrations above 10 ng/mL. Nifoxipam, being highly unstable in serum or during alkaline extraction, was the only compound that could not be detected at relevant trace concentrations in serum. Molecular ions of its published metabolites or degradation products^[4] could not be detected in MS¹.

For each analyte, a linear calibration range (cal_{Low} to cal_{High}) was determined, and calculated concentrations within this range are reported as semi-quantitative result in the automatically generated report. Semi-quantitative results in this preliminary study were found to vary within the range that fits the purpose of providing a good

estimation of the concentration with deviations of ± 10 to $\pm 25\%$ at medium concentrations. In contrast to MS/MS approaches, due to the data dependent acquisition of MSⁿ spectra (including the active exclusion of precursors) only MS¹ scan data is available for quantitation. For coeluting compounds, this leads to a higher influence of peak shape and peak area, explaining the relatively high deviations seen in this study.

Nevertheless, this preliminary data demonstrates that semi-quantitative information can be obtained from screening data using single-point calibration. The Toxtyper software automatically processes full scan data from a routine screening approach, thus no modification to the acquisition method is required. Using customized calibration levels and suitable linear ranges, the accuracy obtained allows the discrimination of therapeutic, sub-therapeutic and potentially toxic serum levels.

As previously indicated, use of internal standards is essential for quantifying serum samples, and confirmatory analysis using a validated quantitative approach is mandatory in forensic case work.

Conclusions

The presented method allows automated identification and semi-quantitative determination of 58 benzodiazepines, including 19 designer benzodiazepines. Limits of detection of the assay allow the detection of sub-therapeutic concentrations or concentrations in the low therapeutic range for the majority of medical benzodiazepines, making the screening applicable for clinical research and forensic analysis. The automated semi-quantitative analysis enables a quick toxicological evaluation of the results and aids in the establishment of subsequent analytical strategies in cases of limited sample availability. Although this approach requires a more time consuming sample preparation when compared to routine immunoassays, unambiguous identification and semi-quantitative determination of compounds offers more detailed sample information.



Learn More

You are looking for further Information?
Check out the link or scan the QR code for more details.

www.bruker.com/customer-insight-boettcher



References

- [1] Demme et al., *Poster presentation*, TIAFT Meeting 2005
- [2] Meyer et al., *Poster presentation (TP29)*, 61th ASMS Conference 2013
- [3] Matuszewski et al., *Anal. Chem.*, **2003**, 75, 3019-3030
- [4] Meyer et al., *Anal Bioanal Chem.*, **2016**, 408:3571-3591

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993

ms.sales.bdal@bruker.com – www.bruker.com