

GPC/SEC-MALLS analysis of Hyaluronic Acid

Application Note Medical Analysis

Authors

Günter Reinhold and Friedhelm Gores
contact: greinhold@pss-polymer.com and fgores@pss-polymer.com

Introduction

Hyaluronic acid is a polymer of disaccharides and can have up to 25,000 disaccharide repeat units. The molar mass range is from 5000 Da to 20 000 000 Da. Hyaluronic acid is mainly used in medical and cosmetic applications.

GPC/SEC is the method of choice for measuring the molar mass distribution of hyaluronic acid. True molar masses can be obtained when on-line multi angle light scattering is used in combination with a refractive index detector. The refractive index increment (dn/dc) for light scattering data evaluation can be determined on-line, offline using dedicated instrumentation or taken from literature.

System Requirements

	Conditions
Pump	PSS SECcurity GPC1260 isocratic pump <ul style="list-style-type: none"> flow rate [mL/min]: 0.5 mobile phase: phosphate buffered saline pH 7.4 (0.005M phosphate buffered saline, 0.069 M NaCl)
Injection system	PSS SECcurity GPC1260 Autosampler <ul style="list-style-type: none"> injection volume 100μL
Columns	<ul style="list-style-type: none"> PSS SUPREMA precolumn (8*50mm) PSS SUPREMA 10μ 30Å (8*300mm) PSS SUPREMA 10μ 10 000Å (8*300mm) PSS SUPREMA 10μ 10 000Å (8*300mm)
Loading	0.5 mg/mL, 100 μ L injection volume
Detectors	Refractive index PSS SECcurity 1260 RI 7 angle MALLS detector 638nm
Software	PSS WinGPC UniChrom optional: Compliance Pack



Procedure, Results & Discussion

High molar mass samples require lower flow-rates and concentrations in GPC/SEC. The GPC/SEC conditions and columns have been optimized for high molar mass hyaluronic acid with respect to loading, flow-rate and column particle size and porosity. For sample preparation the water content of approx. 12 % has been taken into account.

Due to the lack of high molar mass calibration standards the use of light scattering detection is recommended. This allows also to measure true molar masses.

A refractive index increment (dn/dc) for hyaluronic acid of 0.165 has been used to evaluate the

light scattering data.¹ Figure 1 shows the slice concentration measured by the RI detector as well as the on-line measured molar mass. Sample recovery (compare Conc. Calculation vs. Given) was nearly 100% showing that all material eluted from the column.




Figure 2 shows the mass distribution including the cumulative distribution of a hyaluronic acid.




Table 1 shows the molar mass averages of the hyaluronic acid sample obtained using a Pullulan calibration (RI only, apparent molar masses) and using light scattering. The light scattering molar masses are 3 times lower than that obtained with a pullulan calibration curve.

	<i>Pullulan Calibration*</i>	<i>MALLS</i>
Mn [Da]	1 224 000	932 000
Mw [Da]	3 818 000	1 200 000
Mp [Da]	1 398 000	933 000

* Highest molar mass standard: Mp = 2 500 000 Da, P/N PSS-dpul2.5m

Literature:

¹ Lavrenko, Linow, Gornitz in Analytical Ultracentrifugation in Biochemistry and Polymer Science 517-531 Royal Society of Chemistry, Cambridge (1992)