


IC Application Note M-8

Determination of chromate in water applying IC-ICP/MS detection



Hexavalent chromium (chromate) is considered toxic and potentially carcinogenic for which reason its concentration in environmental and drinking water should be as low as possible. The determination of Cr(VI) can be ideally done by ion chromatography coupled to ICP/MS. The separation is performed on a Metrosep A Supp 1 Guard/4.6 column. The chromatogram stems from a chromium-spiked (1 µg/L) wastewater sample.

Results

| | Conc. [µg/L] | SIM ion [m/z] |
|--------------------|-----------------|------------------|
| Waste water spiked | 1.0 | 52 |

Sample

Wastewater spiked with 1 µg/L chromate

Sample preparation

None

Columns

| | |
|-----------------------------|------------|
| Metrosep A Supp 1 Guard/4.6 | 6.1005.340 |
|-----------------------------|------------|

IC Solutions

| | |
|-----------------------|-----------------------------|
| Eluent | 4.6 mmol/L sodium carbonate |
| Suppressor regenerant | 100 mmol/L sulfuric acid |
| Rinsing | Ultrapure water |

ICP/MS Solutions

| | |
|---------------------|-------------------|
| Internal standard 1 | 0.1 mg/L yttrium |
| Internal standard 2 | 0.1 mg/L scandium |

Parameters

| | |
|------------------|-------------|
| Flow rate | 0.35 mL/min |
| Injection volume | 500 µL |
| P _{max} | 15 MPa |
| Recording time | 20 min |

Parameters ICP/MS

| | |
|-------------------|---------------|
| Torque | 2.5 mm |
| Nebulizer | 0.4 mL/min |
| Internal Std Flow | 0.25 mL/min |
| Mode | Collision |
| M/z | 52 |
| Sea Spray | Gacp |
| Spray Chamber | Double pass |
| Acquisition mode | Time resolved |

Analysis

ICP/MS detection

Instrumentation

| | |
|--------------------------------|------------|
| 930 Compact IC Flex ChS/PP/Deg | 2.930.1360 |
| ICP-MS Agilent 7700 | |
| Cetac ASX Autosampler | |
| Remote box | 6.2148.010 |



www.metrohm.com

 **Metrohm**