

What is thermal desorption?

A guide to the history, principles and applications of thermal desorption for gas chromatography

www.markes.com



What is thermal desorption?

“

Thermal desorption is a versatile pre-concentration technique for gas chromatography that is used to analyse volatile and semi-volatile organic compounds.

”



In this e-book, we'll explain:

01 The basics of TD

How it evolved, how it works, why it's so versatile, and aspects of sorbent choice

02 Sampling options

Pumped and passive, canisters and on-line, headspace and more

03 Applications

Why do analysts use TD, and what do they use it for?

We hope that, by the end of this document, you'll have an understanding of how TD benefits GC analysts across a wide range of scenarios – and how it might be able to help you.

For more about TD, visit our website at www.markes.com

WHAT IS THERMAL DESORPTION?

01

THE BASICS



 Return to
contents

Why use TD?

Thermal desorption arose out of the need to improve upon conventional sample preparation techniques for gas chromatography (GC).

GC is very popular for analysing mixtures of organic compounds, but getting the analytes from the sample and into the vapour phase can be difficult.


Thermal desorption overcomes the drawbacks of conventional sample preparation techniques such as solvent extraction, solid-phase micro-extraction, purge-and-trap, and static headspace.

Specifically, it:

- Gives greater sensitivity than these techniques.
- Can be used for a wider range of compound classes.
- Is applicable to a wide range of sample types.
- Is safer and more environmentally-friendly than solvent extraction.
- Is easily automated.
- Is easy to validate.
- Complies with key standard methods.

How does TD work?

The majority of TD applications use sorbent tubes and a two-stage desorption process to focus the analytes into a narrow band of gas, and so achieve the maximum sensitivity enhancement.



The versatility of TD

A major advantage of thermal desorption is that it can be applied to such a wide range of analytes and samples.


RANGE OF ANALYTES

Organic compounds
boiling down to
-100°C...




C₂ hydrocarbons

...to compounds
boiling up to
500°C...



Polycyclic aromatic hydrocarbons

...including reactive
species



Sulfur compounds

RANGE OF SAMPLE TYPES


From



Gases

- Pumped sampling
- Passive sampling
- On-line
- Canisters


...to...



Liquids

- Immersive sorptive extraction

...to...



Solids

- Dynamic headspace
- Headspace sorptive extraction
- Direct desorption

Choosing TD sorbents

Sorbents lie at the heart of why thermal desorption works. Packing tubes and traps with the appropriate sorbent(s) is therefore crucial for the success of all TD methods.

Sorbent type:

The porous polymer sorbent Tenax® TA is the most popular sorbent for TD

Porous polymers



Graphitised carbon blacks



Carbonised molecular sieves






Zeolite molecular sieves



Sorbent strength:

The weaker the sorbent, the better it is at analysing heavier and reactive molecules.



Strong sorbents are usually used for monitoring small, volatile analytes.

Water retention:

Weaker sorbents are generally more hydrophobic, making them a better choice for sampling humid environments.



Conversely, strong sorbents tend to be highly hydrophobic.

The evolution of TD

Over the course of nearly five decades, the capability of thermal desorption has advanced greatly, aided by advances in equipment.

Early 1970s: Scientists begin to experiment with thermal desorption by packing standard GC injector liners with sorbent material.

Mid-1970s: Environmental Monitoring Systems introduces the first commercial thermal desorber.

1981: PerkinElmer introduces the first commercial automated thermal desorber.

Early 1990s: Technical improvements made to improve performance.



1998: Markes introduces the UNITY™ thermal desorber.

2001: Markes introduces ULTRA™ and Air Server™ modules to extend the sampling options.



2005: Markes introduces the TT24-7™ thermal desorber for continuous monitoring.



2008: Markes launches series 2 thermal desorbers.



Late 2000s: Markes introduces range of accessories to expand the application range of TD.



2016: Markes introduces the 'xr' series of thermal desorbers.



1970

1980

1990

2000

2010

2020

3½-inch × ¼-inch TD tubes made standard

Two-stage TD
Electrical cooling
Pre-desorption checks
50-tube automation

Efficient focusing trap
Low-volume valving
Tube sealing for automation

Quantitative re-collection of split flows
Backflush desorption

Analysis of tubes, canisters and on-line air streams
Internal standard addition
100-tube automation

Twin-trap operation for continuous sampling

Analysis of ultra-volatiles, high-boilers and reactive compounds
Automation for tube, canister and on-line monitoring
RFID sample tracking

TD interfaces for dynamic headspace, sorptive extraction, breath sampling, and more

Fully automated re-collection of inlet and outlet split flows
Quantitative analysis of C₂ to C₄₄

WHAT IS THERMAL DESORPTION?

02

SAMPLING OPTIONS




Return to
contents

Sampling options for TD

Thermal desorption can be used with a range of sample introduction methods. Each of those mentioned below is described on the following pages.

Sorbent tubes



Dynamic headspace



Breath



These methods use two stages of sorbent trapping, and so are ideal for trace-level analysis.




Image credits: 1 SP Technical Research Institute of Sweden. 2 Owlstone Medical. 3 Equipco.

Pumped sampling

BASIC OPERATION

- Air is pulled through a sorbent tube using a pump
- Single-bed or multi-bed tubes can be used
- A constant flow of gas into the vessel aids release of analytes onto a sorbent tube

TYPICAL ANALYTES MONITORED

- C₃ to C₄₄

PROS AND CONS

- + Sampling is quick (minutes to hours)
- + A wide range can be monitored in one run
- Sampling conditions need to be optimised so that the most volatile compounds don't 'break through' the sorbent bed


COMMONLY USED FOR

- All types of ambient air monitoring



TOP TIP

A range of tubes are available ready-optimised for key applications



Passive sampling

BASIC OPERATION

- Analytes migrate across an air gap onto a bed of sorbent
- Only single-bed tubes can be used

TYPICAL ANALYTES MONITORED

- C_3 to C_{44}

PROS AND CONS

- + Easy to deploy
- A limited range can be monitored in one run
- Relatively slow sampling


COMMONLY USED FOR

- Factory fenceline monitoring
- Personal monitoring



TOP TIP

Quantitation is easy if you know the uptake rate – and these are available for sampling of a large number of compounds onto common sorbents using industry-standard $3\frac{1}{2}'' \times \frac{1}{4}''$ tubes



Canister sampling

BASIC OPERATION

- Air is allowed to flow at a controlled rate into an evacuated canister
- The air sample is then transferred to a sorbent-packed focusing trap

TYPICAL ANALYTES MONITORED

- Freons
- C₂ hydrocarbons up to C₁₀

PROS AND CONS

- + Ideal for ultra-volatiles
- Not suitable for SVOCs
- Single-stage desorption restricts detection limits
- Bulky to transport


COMMONLY USED FOR

- Ambient air monitoring for volatile 'air toxics'
- Atmospheric research



TOP TIP

In-line water removal devices greatly improve analytical performance



On-line monitoring

BASIC OPERATION

- Air is pulled directly into a sorbent-packed focusing trap

PROS AND CONS

- + Ideal for ultra-volatiles
- Not suitable for SVOCs
- Single-stage desorption restricts detection limits

TYPICAL ANALYTES MONITORED

- Freons
- C₂ hydrocarbons up to C₁₄


COMMONLY USED FOR

- Continuous monitoring for source apportionment
- Rapid detection of hazardous compounds
- Atmospheric research



TOP TIP

Some systems allow continuous monitoring (using two traps working alternately)



Dynamic headspace sampling

BASIC OPERATION

- The sample is placed in a headspace vessel – options include environmental chambers ($\sim 1 \text{ m}^3$), microchambers (44 or 114 cm^3), sampling bags and cuvettes
- A constant flow of gas into the vessel aids release of analytes onto a sorbent tube

TYPICAL ANALYTES MONITORED

- Generally C_3 to C_{20}

PROS AND CONS

- + Highly sensitive
- + Easily adapted to the requirements of the application
- Standard methods for environmental chambers are often time-consuming


COMMONLY USED FOR

- Product certification/compliance
- Quality control
- Rapid product profiling for R&D



TOP TIP

Microchambers can be fitted with accessories for surface and permeation testing



Sorptive extraction

BASIC OPERATION

A polymer-based sorbent is:

- Placed in a liquid sample (*immersive sorptive extraction*)
- Suspended above it (*headspace sorptive extraction*)
- The sorbent is then placed in an empty TD tube and desorbed

TYPICAL ANALYTES MONITORED

- Generally C₃ to C₃₀

PROS AND CONS

- + Suitable solids and liquids
- + Considerably easier than solvent extraction
- Analytes need to be relatively concentrated


COMMONLY USED FOR

- Aroma profiling of foods and beverages



TOP TIP

Immersive sampling is better than headspace sampling for less volatile compounds



Direct desorption



BASIC OPERATION

- A small sample is heated in an empty TD tube

PROS AND CONS

- + Quick and easy
- + Suitable for solids and semi-solids
- Only applicable to relatively homogeneous samples

TYPICAL ANALYTES MONITORED

- Generally C₃ to C₃₀


COMMONLY USED FOR

- Analysis of solids, pastes and resins
- Rapid profiling of polymers
- Analysis of car trim



TOP TIP

Glass tubes make it easy to assess the state of the sample



Breath sampling

BASIC OPERATION

- Breath is captured onto a sorbent tube

PROS AND CONS

- + Less invasive than sampling bodily fluid
- + Sampling doesn't require medical professionals
- Natural variability in breath makes analysis challenging

TYPICAL ANALYTES MONITORED

- Generally C₄ to C₂₀

COMMONLY USED FOR

- Occupational health monitoring
- Research into disease diagnosis



TOP TIP

Devices that capture the last portion of air from the lungs provide the most useful results



 Return to
Sampling Options

WHAT IS THERMAL DESORPTION?

03

APPLICATIONS

 Return to
contents

Why do analysts use TD?

“The need for advanced analytical techniques like thermal desorption is driven by a number of factors.”



Regulations:

Many national and international regulations cite standard methods that rely on TD to achieve the necessary analytical performance.

Emerging products:

Public concern over the health impacts of everyday products – such as fragranced goods – demands rigorous, reliable information on the chemicals they contain. TD is a valuable item in the analyst's toolkit for investigating such products.

Lowering limit levels:

Our understanding of the long-term impacts of chemicals on our health and the environment is constantly improving. Acceptable levels of these chemicals are falling as a result, increasing the need for highly sensitive analytical techniques like TD.

Quality control:

Manufacturers are continually looking for ways to improve their products... for example, by improving consistency, or understanding the causes of off-odours that give rise to customer complaints.

Brand comparison:

The repeatability of thermal desorption and associated sampling techniques make it valuable for carrying out rigorous product comparisons.

The versatility of thermal desorption makes it applicable to a wide range of areas.

For examples of TD in action, download one of our Application Guides.

MARKES
international


Thermal Desorption Applications Guide:

Environmental monitoring


A comprehensive guide to monitoring chemicals in the environment and the workplace using thermal desorption



**Ambient air
monitoring**




**Automotive
studies**



**Soil gas
and water
monitoring**



**Human
health**



**Consumer
environmental
health**




**Food and
drink**




Forensics



**Fragrance
and odour
profiling**




**Industrial air
monitoring and
occupational
health**



**Chemical
ecology**




**Defence and
homeland
security**



**Ecosystems and
the environment**

Our customers

Companies and organisations around the world use Markes' equipment for product quality control and research.



We'd better point out that although we supply all the companies listed on this page, we're not indicating that they've provided an actual endorsement or recommendation. We like to think they might though!

Want to know more?

If you'd like to know more about
the technical aspects of thermal
desorption, or about our products,
you'll find everything you need to know
on our website,

www.markes.com

Alternatively, you'll find our contact details
on the back page – we're always happy
to answer your questions.

MARKES
international

Markes International

Since 1997, Markes International has been at the forefront of innovation for enhancing the measurement of trace-level volatile and semi-volatile organic compounds (VOCs and SVOCs) by gas chromatography (GC).

Our range of thermal desorption products has for many years set the benchmark for quality and reliability. By lowering detection limits, and increasing the options open to the analyst, our thermal desorbers greatly extend the application range of GC.

Our comprehensive portfolio of thermal desorption products includes instruments such as UNITY-xr™ and TD100-xr™, a wide range of high-quality sorbent tubes, and innovative accessories that allow representative vapour profiles to be collected with minimal inconvenience.

Markes is headquartered near Cardiff, UK, and also has laboratory and demonstration facilities in Sacramento, USA, and near Frankfurt, Germany. Markes is a company of the Schauenburg International Group.

Markes International Ltd

Gwaun Elai Medi-Science Campus, Llantrisant, RCT CF72 8XL, UK

T: + 44 (0)1443 230935 **F:** + 44 (0)1443 231531 **E:** enquiries@markes.com

Markes International, Inc.

2355 Gold Meadow Way, Gold River, Sacramento, California 95670, USA

T: 866-483-5684 (toll-free) **F:** 513-745-0741 **E:** enquiries@markes.com

Markes International GmbH

Schleussnerstrasse 42, D-63263 Neu-Isenburg, Frankfurt, Germany

T: +49 (0)6102 8825569 **F:** +49 (0)6102 8825583 **E:** enquiries@markes.com

www.markes.com

