

Method Development for Reversed Phase Chiral LC/MS/MS Analysis of Stereoisomeric Pharmaceutical Compounds with Polysaccharide-based Stationary Phases

Philip J. Koerner, Kari Carlson, Liming Peng, Tivadar Farkas et al.
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Five different Lux® polysaccharide-based chiral stationary phases were explored in the reversed phase elution mode using mobile phases consisting of 0.1 % formic acid in acetonitrile or methanol to demonstrate the feasibility of LC/MS/MS analysis of a variety of acidic pharmaceutical racemates.

Introduction

Developing simple and straightforward reversed phase chiral LC separations coupled with highly sensitive MS detection is a challenging requirement for conducting drug metabolism and pharmacokinetic studies of stereoisomers. Polysaccharide derivatives are the most widely used chiral stationary phases (CSP) due to their wide chiral recognition and high loading capacity. As normal phase is favorable for their principal mechanism of chiral recognition – hydrogen bonding interaction – the majority of chiral separations with polysaccharide phases are performed in normal phase using hexane and alcohol modifiers as mobile phase components. However, these mobile phases are highly flammable and are not compatible with atmospheric pressure ionization (API) MS ion sources. The current research describes the effectiveness of an acidic mobile phase for the separation and detection of using acidic stereoisomers by ESI or APCI LC/MS/MS and for method development of these applications.

Experimental

Analytes: 15 acidic compounds of pharmaceutical interest were analyzed. The structures are shown in **Figure 2** and the MS ionization and mass reaction monitoring (MRM) transitions monitored are listed in **Table 1**.

- Columns:** Lux 3 μ m Cellulose-1, 150 x 2.0 mm
 Lux 3 μ m Cellulose-2, 150 x 2.0 mm
 Lux 3 μ m Amylose-2, 150 x 2.0 mm
 Lux 5 μ m Cellulose-1, 250 x 4.6 mm
 Lux 5 μ m Cellulose-4, 250 x 4.6 mm
 Lux 5 μ m Cellulose-3, 250 x 4.6 mm
 Lux 5 μ m Amylose-2, 250 x 4.6 mm
 Kinetex® 2.6 μ m C18, 50 x 2.1 mm (used for achiral analysis)
- Flow Rate:** 0.2 mL/min (3 μ m, 150 x 2.0 mm) or
 1.0 mL/min (5 μ m, 250 x 4.6 mm) – flow split to 0.25 mL/min into MS/MS
- Temperature:** 25 °C
- Detection:** UV @ 220 nm
- Injection Volume:** 5 μ L (150 x 2.0 mm) or 10–20 μ L (250 x 4.6 mm)
- Mobile Phases:** 1. 0.1 % Formic acid in Acetonitrile or Methanol
 2. 5 mM Ammonium bicarbonate in Acetonitrile or Methanol (achiral analysis)
 3. 5 mM Ammonium formate in Acetonitrile or Methanol (achiral analysis)
 4. 5 mM Ammonium acetate in Acetonitrile or Methanol (achiral analysis)
- Instrument:** HPLC System: Agilent® 1200 series equipped with binary pump and autosampler (Agilent, Palo Alto, CA)
 MS Detector: AB SCIEX™ 4000 LC/MS/MS
 Turbo V™ source with ESI or APCI probe
 MS Detection: TurbolonSpray® – ESI or APCI in Positive or Negative Ion Mode; MRM

Results and Discussion

Chiral LC/MS/MS Experiments

Five different polysaccharide-based chiral stationary phases (Lux Cellulose-1, Lux Cellulose-2, Lux Cellulose-3, Lux Cellulose-4, and Lux Amylose-2 (**Figure 1**)) were explored in the reversed phase elution mode for the separation of a variety of acidic compounds of pharmaceutical interest in mobile phases consisting of 0.1% formic acid in acetonitrile or methanol and MS/MS detection.

Figure 1.
Structures of Polysaccharide-based Chiral Stationary (CSPs) Phases





Table 1.
MRM Transitions and Concentrations of Acidic Racemates

Compound	IS and MRM	Conc.*	Compounds	IS and MRM	Conc.*
Ibuprofen	ESI 205.1/160.1	100	Abscisic acid	APCI 263.0/152.7	50
Flurbiprofen	ESI 243.0/198.7	50	Mecoprop	ESI 214.1/141.7	50
Suprofen	ESI 261.1/111.0	50	Ketorolac	ESI 254.0/209.8	50
Fenoprofen	APCI 241.0/196.8	100		ESI 256.2/105.0	
Carprofen	APCI 272.8/228.8	100	Etodolac	ESI 286.1/242.0	50
Indoprofen	ESI 282.1/236.1	50	Warfarin	ESI 307.1/160.9	50
Proglumide	ESI 333.1/120.9	50		ESI 309.2/163.1	
1-(Phenylsulfonyl)-3-indoleboronic acid	ESI 300.0/209.8	100	Bendroflumethiazide	ESI 420.1/77.9	50
			Trichlormethiazide	ESI 377.8/241.6	50

* Conc. (ng/mL)

Figure 2.
Molecular Structures of Acidic Racemates



Selection of Mobile Phase Additives

As acidic analytes are present as anions in mobile phases of neutral pH their retention is not favorable on polysaccharide-based CSPs. Early elution and poor or no enantioseparation can result under these conditions. Acidic mobile phase additives are often required to suppress the dissociation of acidic analytes, resulting in increased retention and improved enantioselectivity.

Three volatile organic acids, (trifluoroacetic acid (TFA, pK_a 0.59), formic acid (FA, pK_a 3.75), and acetic acid (HAc, pK_a 4.76) were evaluated on Lux® Cellulose-1 and Lux Amylose-2 CSPs as acidic additives. In general, these additives provide similar enantioresolution for weakly acidic racemates (**Figures 3A-3C**) while the stronger acidic additive (TFA) performs better for the stronger acidic racemates (**Figures 3D-3F**). Formic acid provides comparable enantioseparations and peak shapes to TFA. Considering the “memory effect” of TFA commonly experienced on polysaccharide-based CSPs and its ion suppressing tendency in MS detection, formic acid was selected in preference over TFA in reversed phase mobile phase for the chiral separation and MS/MS detection of acidic racemates.

Effect of Mobile Phase Additives

The LC/MS/MS responses of acidic racemates with ESI negative (ESI-) mode in 5 mM ammonium formate, 5 mM ammonium acetate, and 5 mM ammonium bicarbonate containing mobile phases with either acetonitrile or methanol as organic modifier using an achiral column (Kinetex® 2.6 μ m C18) were compared to responses in 0.1% formic acid (**Figures 4 and 5**). The results show that MS/MS responses using 0.1% formic acid in acetonitrile or methanol are comparable to the responses using the other acid additives in ESI- mode, with the exception of carprofen which showed poor response with all of the acidic mobile phase additives. This shows that using 0.1% formic acid as acidic mobile phase additive is fully compatible with MS/MS detection and can be implemented as the first choice for mobile phase additive.

Effect of Organic Modifier on Chiral Resolution

Acetonitrile or methanol was the organic modifier used in chiral reversed phase HPLC. Decreasing the eluting strength of the mobile phase by decreasing the percentage of acetonitrile or methanol in the mobile phase will increase retention and resolution (**Figure 6**). However, once enantiomers elute later than about 10 minutes with only partial resolution, baseline separation can be rarely achieved by further decreasing the % organic modifier in the mobile phase. In our study, acetonitrile was more successful than methanol in providing chiral resolution on Lux CSPs in reversed phase (RP) mode (**Figure 8**), typically yielding sharper and narrower peaks.

TN-1079

APPLICATIONS

Figure 3.
Effect of Acidic Additives on the Enantioseparation of Acidic Racemates in Reversed Phase





Figure 4.
LC/MS/MS Responses of Acidic Racemates in ESI- with Different Additives and Methanol



TN-1079

APPLICATIONS

Figure 5.
LC/MS/MS Responses of Acidic Racemates in ESI- with Different Additives and Methanol





Figure 6.
Effect of Acidic Additives on the Enantioseparation of Acidic Racemates in Reversed Phase

Proglumide

Lux® 3 µm Cellulose-1
Flow Rate: 0.2 mL/min
Dimensions: 150 x 2.0 mm

Ibuprofen

Lux 5 µm Cellulose-3
Flow Rate: 1 mL/min
Dimensions: 250 x 4.6 mm



Column: As noted
Dimensions: As noted
Flow Rate: As noted
Temperature: 25 °C
Detection: Mass Spectrometer (MS)
Injection Volume: 5-20 µL
Mobile Phase: As noted
Sample: As noted

TN-1079


APPLICATIONS

Chiral LC/MS/MS Applications

Figures 7-11 demonstrate 15 chiral separations on several different Lux chiral stationary phases. APCI negative (APCI-) mode was employed for the MS/MS detection of carprofen, abscisic acid, and fenoprofen as it provides much better MS/MS signals than ESI- for these acidic racemates. ESI positive (ESI+) mode was used for the MS/MS detection of suprofen and gave slightly im-


proved signal intensity for warfarin and ketorolac versus ESI- (**Figure 10**). Most compounds evaluated here eluted in less than 10 min with baseline resolution in mobile phases of various elution strength. The results show that Lux® Cellulose-3 was most successful in separating acidic racemates (ten racemates), especially nonsteroidal anti-inflammatory drugs.

Figure 7.
Enantioseparations in Reversed Phase on Lux 3 µm Cellulose-1 with Acetonitrile



Column: Lux 3 µm Cellulose-1
Dimensions: 150 x 2.0 mm
Part No.: 00F-4458-B0
Flow Rate: 0.2 mL/min
Temperature: 25 °C
Detection: Mass Spectrometer (MS)
Injection Volume: 5 µL
Mobile Phase: As noted
Sample: As noted

Figure 8.
Enantioseparations on Lux 3 µm Cellulose-2 with Acetonitrile or Methanol as Modifier



Column: Lux 3 µm Cellulose-2
Dimensions: 150 x 2.0 mm
Part No.: 00F-4458-B0
Flow Rate: 0.2 mL/min
Temperature: 25 °C
Detection: Mass Spectrometer (MS)
Injection Volume: 5 µL
Mobile Phase: As noted
Sample: As noted

TN-1079

APPLICATIONS

Figure 9.
Enantioseparations on Lux® 5 µm Cellulose-4 with Acetonitrile






Figure 10.
Enantioseparations in Reversed Phase on Lux 3 µm Amylose-2 with Acetonitrile



TN-1079

APPLICATIONS

Figure 11.
Enantioseparations in Reversed Phase on Lux® 5 µm Cellulose-3 with Acetonitrile or Methanol as Modifier



Conclusions

The chiral LC/MS/MS analysis of fifteen different acidic racemates are successfully demonstrated on the polysaccharide-based CSPs Lux Cellulose-1, Lux Cellulose-2, Lux Cellulose-3, Lux Cellulose-4, and Lux Amylose-2 in reversed phase (RP) elution mode.

Formic acid is a good first choice for an acidic RP mobile phase additive as it leads to increased retention and improved enantioselectivity for acidic enantiomers and is also compatible with MS/MS detection.

Increasing the percentage of organic modifier (acetonitrile or methanol) in the RP mobile phase has the expected effect of decreasing retention and enantioselectivity. Adjusting the organic modifier content of the mobile phase is therefore essential to optimizing chiral resolution.

Free Chiral Screening Services



Visit us online at: www.phenomenex.com/Phenologix

TN-1079

APPLICATIONS

Australia

t: 02-9428-6444
f: 02-9428-6445
auinfo@phenomenex.com

Austria

t: 01-319-1301
f: 01-319-1300
anfrage@phenomenex.com

Belgium

t: +31 (0)30-2418700
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: (800) 543-3681
f: (310) 328-7768
info@phenomenex.com

Denmark

t: 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: 09 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: 01 30 09 21 10
f: 01 30 09 21 11
franceinfo@phenomenex.com

Germany

t: 06021-58830-0
f: 06021-58830-11
anfrage@phenomenex.com

India

t: 040-3012 2400
f: 040-3012 2411
indiainfo@phenomenex.com

Ireland

t: 01 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: 051 6327511
f: 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: 030-2418700
f: 030-2383749
nlinfo@phenomenex.com

New Zealand

t: 09-4780951
f: 09-4780952
nzinfo@phenomenex.com

Norway

t: 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
f: (310) 328-7768
info@phenomenex.com

Sweden

t: 08 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

Lux® Ordering Information

3 µm Analytical Columns (mm)

	50 x 2.0	150 x 2.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	SecurityGuard™ Cartridges (mm)	
Phases							4 x 2.0*	4 x 3.0*
Cellulose-4	00B-4490-B0	00F-4490-B0	00B-4490-E0	00D-4490-E0	00F-4490-E0	00G-4490-E0	AJ0-8626	AJ0-8627
Cellulose-3	00B-4492-B0	00F-4492-B0	00B-4492-E0	00D-4492-E0	00F-4492-E0	00G-4492-E0	AJ0-8621	AJ0-8622
Cellulose-2	00B-4456-B0	00F-4456-B0	00B-4456-E0	00D-4456-E0	00F-4456-E0	00G-4456-E0	AJ0-8398	AJ0-8366
Cellulose-1	00B-4458-B0	00F-4458-B0	00B-4458-E0	00D-4458-E0	00F-4458-E0	00G-4458-E0	AJ0-8402	AJ0-8403
Amylose-2	00B-4471-B0	00F-4471-B0	00B-4471-E0	00D-4471-E0	00F-4471-E0	00G-4471-E0	AJ0-8471	AJ0-8470

for ID: 2.0–3.0 mm 3.2–8.0 mm

5 µm Analytical Columns (mm)

	50 x 2.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	SecurityGuard™ Cartridges (mm)	
Phases						4 x 2.0*	4 x 3.0*
Cellulose-4	00B-4491-B0	00B-4491-E0	00D-4491-E0	00F-4491-E0	00G-4491-E0	AJ0-8626	AJ0-8627
Cellulose-3	00B-4493-B0	00B-4493-E0	00D-4493-E0	00F-4493-E0	00G-4493-E0	AJ0-8621	AJ0-8622
Cellulose-2	00B-4457-B0	00B-4457-E0	00D-4457-E0	00F-4457-E0	00G-4457-E0	AJ0-8398	AJ0-8366
Cellulose-1	00B-4459-B0	00B-4459-E0	00D-4459-E0	00F-4459-E0	00G-4459-E0	AJ0-8402	AJ0-8403
Amylose-2	00B-4472-B0	00B-4472-E0	00D-4472-E0	00F-4472-E0	00G-4472-E0	AJ0-8471	AJ0-8470

for ID: 2.0–3.0 mm 3.2–8.0 mm

5 µm Semi-Prep Columns (mm)

	150 x 10.0	250 x 10.0	SecurityGuard™ Cartridges (mm)	
Phases			10 x 10.0*	/3pk
Cellulose-4	00F-4491-N0	00G-4491-N0	AJ0-8628	
Cellulose-3	00F-4493-N0	00G-4493-N0	AJ0-8623	
Cellulose-2	00F-4457-N0	00G-4457-N0	AJ0-8399	
Cellulose-1	00F-4459-N0	00G-4459-N0	AJ0-8404	
Amylose-2	00F-4472-N0	00G-4472-N0	AJ0-8472	

for ID: 9–16 mm

5 µm Axia™ Packed Preparative Columns (mm)

	150 x 21.2	250 x 21.2	250 x 30	250 x 50	SecurityGuard™ Cartridges (mm)	
Phases					15 x 21.2**	15 x 30.0*
Cellulose-4	00F-4491-P0-AX	00G-4491-P0-AX	00G-4491-U0-AX	00G-4491-V0-AX	AJ0-8629	AJ0-8630
Cellulose-3	00F-4493-P0-AX	00G-4493-P0-AX	00G-4493-U0-AX	00G-4493-V0-AX	AJ0-8624	AJ0-8625
Cellulose-2	00F-4457-P0-AX	00G-4457-P0-AX	00G-4457-U0-AX	00G-4457-V0-AX	AJ0-8400	AJ0-8401
Cellulose-1	00F-4459-P0-AX	00G-4459-P0-AX	00G-4459-U0-AX	00G-4459-V0-AX	AJ0-8405	AJ0-8406
Amylose-2	00F-4472-P0-AX	00G-4472-P0-AX	00G-4472-U0-AX	00G-4472-V0-AX	AJ0-8473	AJ0-8474

for ID: 18–29 mm 30–49 mm

*SecurityGuard Analytical Cartridges require holder, Part No. : KJ0-4282

*SemiPrep SecurityGuard Cartridges require holder, Part No. : AJ0-7220

**PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8223

*PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8277



If Lux analytical columns (≤ 4.6 mm ID) do not provide at least an equivalent or better separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Lux and Kinetex are registered trademarks, and Axia and SecurityGuard are trademarks of Phenomenex. TurbolonSpray is a registered trademark and AB SCIEX, API 3000, and Turbo V are trademarks of AB Sciex Pte. Ltd.

Disclaimer

Comparative separations may not be representative of all applications.

© 2012 Phenomenex, Inc. All rights reserved.

United Kingdom

t: 01625-501367
f: 01625-501796
ukinfo@phenomenex.com

United States

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

All other countries: Corporate Office USA

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com