



Improved Recoveries and Lower Background for the Analysis of PAHs in Olive Oil

Using Supelclean™ EZ-POP NP and an SLB®-5ms GC Column

A recent survey implemented by Supelco revealed that 58% of analytical chemists in the food industry find the analysis of compounds in fatty matrices to be highly problematic. Lipid interferences can produce elevated detection limits, contamination of LC and GC systems, and ultimately decrease instrument and column lifetime.

Current extraction and cleanup techniques for fatty samples may be time consuming, expensive, and/or provide insufficient background removal prior to analysis. For this reason, a new approach has been developed. The Supelclean EZ-POP NP, a dual-layer SPE cartridge containing Florisil® and Z-Sep/C18, was designed for the extraction of nonpolar analytes from oil matrices. By means of Lewis acid/base and hydrophobic interactions, fatty matrix interferences are preferentially retained by the cartridge while analytes of interest are eluted. In this way, lipid interferences are removed from the sample.

Experimental

The Supelclean EZ-POP NP and two competitor silica gel SPE cartridges were compared in terms of matrix removal and analyte recovery for the extraction of select polycyclic aromatic hydrocarbons (PAHs) from olive oil. Multiple replicates of both unspiked and spiked (20 ng/g with PAHs) oil samples were processed for each SPE cleanup technique. The cleanup procedures are summarized in **Tables 1 and 2**.

Table 1. Cleanup Procedure for Olive Oil using Supelclean EZ-POP NP

1. Add 10 mL of acetone to the Supelclean EZ-POP NP cartridge (54341-U), and allow it to elute with gravity.
2. Dry the cartridge at -10°C Hg for 10 min.
3. Accurately weigh 0.500 mL of oil directly onto the SPE cartridge.
4. Add the internal standard directly to the oil on the cartridge.
5. Allow the sample to penetrate the upper frit.
6. Pulling vacuum, elute the analytes of interest with 15 mL of acetonitrile (rate of approximately 1 drop/sec.).
7. Concentrate the samples to a final volume of 1 mL under nitrogen (5 psi) at 40°C .

Table 2. Cleanup Procedure for Olive Oil using Silica Gel SPE Cartridges (Competitor A and B)*

1. Condition the SPE cartridges with 20 mL of hexane.
2. Dilute 5 g of oil with hexane to a final volume of 10 mL. Load 1 mL onto the cartridge. Add the internal standard directly onto the cartridge.
3. Wash the cartridge with 8 mL hexane:methylene chloride (70:30).
4. Elute the analytes of interest with 8 mL hexane:methylene chloride (70:30).
5. Concentrate the samples to a final volume of 1 mL under nitrogen (5 psi) at 40°C .

* Gravity elution was used for the above steps.

Extracts were then analyzed by GC-MS (SIM mode). Quantitation was performed against a 5-point calibration curve (1–20 ng/mL) prepared in unspiked olive oil extract with naphthalene- d_8 , fluoranthene- d_{10} , perylene- d_{12} internal standards, spiked at 10 ng/mL.

Results and Discussion


Matrix Removal

Figure 1 shows the GC-MS (full scan) chromatograms of the extracts after cleanup with Supelclean EZ-POP NP, Competitor A silica gel SPE, and Competitor B silica gel SPE, respectively. The chromatograms illustrate that the EZ-POP NP cleanup produces much lower background than the silica gel SPE cleanup.


Figure 1. GC-MS Full Scan Chromatograms of Olive Oil Extract (same y axis)

column: SLB®-5ms, 20 m x 0.18 mm I.D., 0.18 μm (28564-U)
oven: 60°C (1 min.), $15^{\circ}\text{C}/\text{min.}$ to 250°C , $8^{\circ}\text{C}/\text{min.}$ to 330°C (7 min.)
inj. temp.: 300°C
carrier gas: helium, 1 mL/min constant flow
detector: MS
MSD interface: 330°C
injection: 1 μL , pulsed splitless (50 psi until 0.75 min, splitter open at 0.75 min.)
liner: 4 mm ID FocusLiner™ with taper and quartz wool

Supelclean EZ-POP NP Cleanup



Competitor A Silica Gel SPE Cleanup



Competitor B Silica Gel SPE Cleanup






Figure 2. Analyte Recovery of PAHs from Olive Oil Extract (n=3)


Analyte Detection and Recovery

As illustrated in **Figure 2**, the EZ-POP NP samples produced better overall analyte recoveries than those cleaned with the silica SPE cartridges. The reproducibility for all PAHs tested, exhibited by percent relative standard deviation (%RSD), fell within an acceptable range of less than 20% for the EZ-POP NP cleaned samples. The EZ-POP NP cleanup demonstrated improvements over silica gel cleanup in the following areas:

- Recoveries for the lighter PAHs (2–3 rings) in the EZ-POP NP cleanup were not affected by matrix as they were for the silica gel samples.
- The heavier PAHs of 4–6 rings showed improved recoveries using EZ-POP NP cartridge over the silica gel cartridges.
- The heavy background from the silica cleanup caused retention time shifts for the last six PAH peaks, making identification difficult. Also, the high level of late eluting background in the Competitor B silica sample prevented accurate quantitation of the last internal standard, perylene- d_{12} . On the contrary, there were no background complications associated with the EZ-POP NP cleanup.

Hence, the EZ-POP NP produced better overall recoveries and removed more problematic matrix interferences than the silica gel SPE cartridges, while maintaining good reproducibility.

Conclusion

This study compared the Supelclean EZ-POP NP to silica gel SPE for the extraction of PAHs from olive oil, a problematic fatty matrix. In terms of lipid removal, the EZ-POP NP removed more unwanted background than silica gel SPE, greatly decreasing the matrix effects. The EZ-POP NP produced better overall analyte recoveries than the silica gel SPE, with adequate reproducibility. Thus, the Supelclean™ EZ-POP NP provides suitable matrix removal for rugged GC-MS analysis of PAHs in olive oil.

Reference

- Moret, S., Conte, L., *J. Sep. Sci.* **25**, 96-100, (2002).

Featured and Related Products

Description	Qty.	Cat. No.
Supelclean EZ-POP NP SPE Cartridge	20	54341-U
Column		
SLB®-5ms Capillary GC Column 20 m x 0.18 mm I.D., 0.18 µm	1	28564-U
Analytical Standards		
Polynuclear Aromatic Hydrocarbons Mix, <i>TraceCERT</i> ® CRM 2000 µg/mL each component in methylene chloride: benzene (1:1)		CRM48905
Naphthalene- d_8 solution 2000 µg/mL in methylene chloride		48715-U
Fluoranthene- d_{10} 50 mg ampul		442843
Perylene- d_{12} solution 2000 µg/mL in methylene chloride		48081
Solvent		
Acetonitrile, LC-MS CHROMASOLV®		34967

Visit our Food & Beverage Analysis web portal at
sigma-aldrich.com/food

Including products from **Fluka**
 Analytical

Enabling Science to
 Improve the Quality of Life

Order/Customer Service: sigma-aldrich.com/order
 Technical Service: sigma-aldrich.com/techservice
 Development/Custom Manufacturing Inquiries **SAFC**® safcglob@ial.com
 Safety-related Information: sigma-aldrich.com/safetycenter

World Headquarters
 3050 Spruce St.
 St. Louis, MO 63103
 (314) 771-5765
sigma-aldrich.com