


# Extraction of Thyroid Hormones, T3, rT3 and T4 from Serum Using EVOLUTE® EXPRESS AX Plates Prior to LC-MS/MS Analysis



**Figure 1.** Structures of T3, rT3 and T4

## Introduction

This application note describes a polymer based-based strong anion exchange mixed-mode SPE protocol for the extraction of T3, rT3 and T4 from serum prior to LC-MS/MS analysis.

The method described in this application note demonstrates selective extraction of the thyroid hormones T3, rT3 and T4 from serum. 500  $\mu$ L of serum was extracted using the EVOLUTE EXPRESS AX 30mg fixed well plate format. High reproducible recoveries and extremely low phospholipid content were observed, demonstrating limits of quantitation of 50 pg/mL.

## Analytes

Tri-iodothyronine (T3), reverse Tri-iodothyronine (rT3) and Thyroxine (T4).

## Sample Preparation Procedure

<b>Format:</b>	<b>EVOLUTE® EXPRESS AX 30 mg Fixed Well Plate, part number 603-0030-PX01</b>
<b>Serum Pre-treatment:</b>	Take 500 $\mu$ L of serum and spike with internal standard (10 $\mu$ L of a 25 ng/mL solution). Add 500 $\mu$ L of a mixture of 5% $\text{NH}_4\text{OH}$ (aq) and acetonitrile (MeCN) (50/50, (v/v)) and vortex mix thoroughly.

## Solid Phase Extraction

<b>Plate Conditioning:</b>	Condition each well with methanol (1 mL)
<b>Plate Equilibration:</b>	Equilibrate each well with water (1 mL)
<b>Sample Loading:</b>	Load pre-treated sample (1 mL) to the 96 well-plate
<b>Wash 1:</b>	Elute interferences with 5% $\text{NH}_4\text{OH}$ (aq)/MeCN (50/50, (v/v), 1 mL)
<b>Wash 2:</b>	Elute interferences with 5% $\text{NH}_4\text{OH}$ (aq) (1 mL)
<b>Wash 3:</b>	Elute interferences with Methanol (1 mL)
<b>Wash 4:</b>	Elute interferences with 2% formic acid in dichloromethane (1 mL)
<b>Elution:</b>	Elute analytes with 5% formic acid in methanol (500 $\mu$ L) into a collection plate containing 2 $\mu$ L of ethylene glycol in each well
<b>Post Elution:</b>	Evaporate at 40 °C in a stream of air or nitrogen using a SPE Dry. <b>Note:</b> The use of ethylene glycol helps to minimize non-specific binding to the plastic collection plate by avoiding complete drying during the evaporation procedure.
<b>Reconstitution:</b>	Reconstitute in 50/50 (v/v) UPLC Mobile Phase A/B (100 $\mu$ L)

## UPLC Conditions

<b>Instrument:</b>	Waters ACQUITY UPLC
<b>Column:</b>	Phenomenex Kinetex C18 UHPLC column (1.7 $\mu$ , 100 x 2.1 mm id)
<b>Mobile Phase:</b>	A: 2 mM ammonium acetate/0.1% formic acid (aq) B: 2 mM ammonium acetate/0.1% formic acid/MeOH
<b>Flow Rate:</b>	0.3 mL/min

**Table 1.** Gradient Conditions

Time	% A	% B	Curve
0	50	50	1
2.5	10	90	6
2.6	50	50	11

Curve 11: Conditions in line initiated immediately once time passed. i.e. 50:50 resumed at 2.6 minutes

Curve 6: Linear Gradient

<b>Injection Volume:</b>	15 $\mu$ L (partial loop with overfill)
<b>Sample Temperature:</b>	20 °C
<b>Columns Temperature:</b>	40 °C

## Mass Spectrometry Conditions

Ions were selected in order to achieve maximum sensitivity, and the MS operated in positive ion polarity mode, using multiple reaction monitoring.

<b>Instrument:</b>	Premier XE triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis
<b>Desolvation Temperature:</b>	450 °C
<b>Ion Source Temperature:</b>	150 °C
<b>Collision Cell Pressure:</b>	3.76 e <sup>-3</sup> mbar


Positive ions acquired in the multiple reaction monitoring (MRM) mode:

**Table 2.** MRM Conditions (Qualifier ion details shown in parenthesis)

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
T3	651.8 > 605.8 (651.8 > 479.0)	35 (35)	21 (35)
rT3	651.8 > 605.8 (651.8 > 508.0)	35 (35)	21 (20)
T3/rT3-d6 ISTD	657.9 > 611.9	35	21
T4	777.8 > 731.8 (777.8 > 351.2)	40 (40)	25 (45)

## Results


Good chromatographic separation of T3/rT3 and T4 was achieved in less than 2.5 minutes, as shown in **Figure 2**.



**Figure 2.** Chromatographic separation of T3, rT3 and T4 from spiked serum at 100 pg/mL

## Recovery

Stripped serum was spiked at various levels from 50-5000 pg/mL. High reproducible recoveries > 80% with corresponding RSDs < 10% were demonstrated. Typical recovery data is shown in **Figure 3**.




**Figure 3.** Recovery profile for Thyroid hormones extracted at 2 ng/mL

## Calibration Curves


Calibration curves were generated using stripped serum spiked at concentrations from 50-5000 pg/mL. Good coefficients of determination were obtained for T3 and rT3 ( $r^2 > 0.99$ ). Endogenous levels of T4 contributing towards a substantial intercept affected calibration curve performance, demonstrating  $r^2$  values of only 0.97.

Compound name: T3 Quant  
Correlation coefficient:  $r = 0.996739$ ,  $r^2 = 0.993489$   
Calibration curve:  $1.138 \times 10^3 x + 9.73202$   
Response type: Internal Std (Ref 7), Area \* (IS Conc / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




**Figure 4.** Spiked serum calibration curve demonstrating T3 linearity from 50-5000 pg/mL (n=4).

Compound name: rT3 Quant  
Correlation coefficient:  $r = 0.996714$ ,  $r^2 = 0.993439$   
Calibration curve:  $0.95792 \times 10^3 x + 0.223514$   
Response type: Internal Std (Ref 8), Area \* (IS Conc / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



**Figure 5.** Spiked serum calibration curve demonstrating rT3 linearity from 50-5000 pg/mL (n=4).


Compound name: T4 Quant  
Correlation coefficient:  $r = 0.95557$ ,  $r^2 = 0.971913$   
Calibration curve:  $1.06787 \times 10^3 x + 746.617$   
Response type: Internal Std (Ref 9), Area \* (IS Conc / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



**Figure 6.** Calibration Curve demonstrating T4 linearity from 50-5000 pg/mL.

## Phospholipid Removal

The robust washing protocols that are possible when using EVOLUTE® EXPRESS AX plates led to extremely clean extracts with very low matrix effects. Residual phospholipids were investigated to provide an indication of extract cleanliness. We investigated the most abundant phospholipids (selected from full scan, SIR and precursor ion scanning experiments) using MRM transitions monitoring the common 184 product ion. Figure 8 demonstrates phospholipid content comparing protein precipitated serum, solvent blank and the final extraction protocol.



**Figure 7.** MRM TICs comparing phospholipid content

## Additional Notes

### 1. Buffer Preparation

- 5%  $\text{NH}_4\text{OH}$  (aq): Take 95 mL of  $\text{H}_2\text{O}$  and add 5 mL of  $\text{NH}_4\text{OH}$  (28–30% concentration).
- 5%  $\text{NH}_4\text{OH}$  (aq)/MeCN (50/50 (v/v)): Take 50 mL of 5%  $\text{NH}_4\text{OH}$  (aq) and add 50 mL of MeCN.
- 2 % formic acid in dichloromethane: Take 98 mL of DCM and add 2 mL of formic acid (98% concentration).
- 5 % formic acid in methanol: Take 95 mL of MeOH and add 5 mL of formic acid (98% concentration).
- 2 mM ammonium acetate/0.1% formic acid (aq): Weigh 0.15416 g and dissolve in  $\text{H}_2\text{O}$ . Add 1 mL of formic acid and make up to 1 L in  $\text{H}_2\text{O}$ .
- 2 mM ammonium acetate/0.1% formic acid/MeOH: Weigh 0.15416 g and dissolve in MeOH. Add 1 mL of formic acid and make up to 1 L in MeOH.

### 2. Processing Conditions

- Positive Pressure: Process at approximately 2–3 psi.
- Vacuum Processing: Process at approximately -0.2 bar.

### 3. Conditioning and equilibration steps

- EVOLUTE EXPRESS 96-well plates may be used without conditioning and equilibration steps. Evaluate the requirement for these steps during method development.

## Ordering Information

Part Number	Description	Quantity
603-0030-PX01	EVOLUTE® EXPRESS AX 30 mg Fixed Well Plate	1
121-9600	Biotage® VacMaster™-96 Sample Processing Manifold	1
PPM-96	Biotage® PRESSURE+ 96 Positive Pressure Manifold 96 Well	1
SD-9600-DHS-EU	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
SD-9600-DHS-NA	Biotage® SPE Dry Sample Concentrator System 100/120 V	1

For the latest application notes visit [www.biotage.com](http://www.biotage.com)

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
[order@biotage.com](mailto:order@biotage.com)  
 Support Tel: +46 18 56 59 11  
 Support Fax: +46 18 56 57 11  
[eu-1-pointsupport@biotage.com](mailto:eu-1-pointsupport@biotage.com)

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
[ordermailbox@biotage.com](mailto:ordermailbox@biotage.com)  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
[us-1-pointsupport@biotage.com](mailto:us-1-pointsupport@biotage.com)

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
[jp\\_order@biotage.com](mailto:jp_order@biotage.com)  
[jp-1-pointsupport@biotage.com](mailto:jp-1-pointsupport@biotage.com)

### CHINA

Tel: +86 21 2898 6655  
 Fax: +86 21 2898 6153  
[cn\\_order@biotage.com](mailto:cn_order@biotage.com)  
[cn-1-pointsupport@biotage.com](mailto:cn-1-pointsupport@biotage.com)

To locate a distributor,  
 please visit our website at  
[www.biotage.com](http://www.biotage.com)

#### Part Number: AN817

© 2014 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. Product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

For more information visit [www.biotage.com](http://www.biotage.com).