

# Analysis of Vitamin D in Milk and Infant Formula using UHPLC/MS/MS per AOAC Method 2011.11

Olga Shimelis, Principal R&D Scientist, Supelco  
olga.shimelis@sial.com

Fortification of milk, dairy products, infant and adult nutritionals with vitamin D is required in many countries to provide additional dietary vitamin D. The fortification of milk in the USA, for example, is done at 400 International Units (IU) per eight ounces of milk. The recently adopted AOAC method 2011.11 describes the procedure for the analysis of vitamin D in fortified infant formula. This method was used in the current work to measure vitamin D in the fortified whole milk and in an infant formula sample. Two forms of vitamin D are recognized: vitamin D<sub>2</sub> or ergocalciferol and vitamin D<sub>3</sub> or cholecalciferol. Although both forms may be present in the fortified food, vitamin D<sub>3</sub> is more commonly used for fortification of dairy products. Vitamin D hydroxy-metabolites can also be present in food products at low levels. While quantitative analysis was performed for vitamins D<sub>2</sub> and D<sub>3</sub>, the samples were scanned for the presence of 25-hydroxy metabolites.

## Method

AOAC method 2011.11 was followed for analysis of vitamin D in both milk and infant formula. Vitamin D<sub>2</sub> and D<sub>3</sub> standards and their deuterated analogues were used. The UHPLC column was Titan C18, 10 cm x 2.1 mm packed with 1.9 µm particles. The HPLC method was adopted from the standard AOAC method, including mobile phases, gradient, and flow rates. Baseline separation between vitamin D<sub>3</sub> and D<sub>2</sub> peaks, and between all vitamins and matrix component peaks was obtained on the Titan C18 UHPLC column. Analysis of vitamin D was done using MS/MS with APCI detection. MS transitions used in this study are shown in Table 1.

**Table 1. MS Parameters Used for the Measurement of Vitamin D Forms**

Compound	Q1	Q3
Vitamin D <sub>3</sub>	385.5	259.2
Vitamin D <sub>3</sub>	385.5	159.2
Vitamin D <sub>2</sub>	397.5	125.2
Vitamin D <sub>2</sub>	397.5	107.2
Vitamin D <sub>3</sub> -d <sub>3</sub>	388.5	259.2
Vitamin D <sub>2</sub> -d <sub>3</sub>	400.5	125.2
25-OH-Vitamin D <sub>3</sub>	383.4	257.2
25-OH-Vitamin D <sub>2</sub>	395.2	269.2

The sample preparation for both milk and infant formula samples was done through saponification, extraction, evaporation, and reconstitution of extracts into LC-compatible solvents. The sample preparation steps following AOAC method 2011.11 are shown in Figure 1.

**Figure 1. Preparation Method for Milk and Nutritional Supplements**

<b>Saponification</b>	<ul style="list-style-type: none"> <li>To 30 mL of milk or infant formula add 40 mL of reagent-grade alcohol containing 2% pyrogallol, isotope D<sub>3</sub> I.S., and 20 mL of 50% w/v KOH</li> <li>Mix overnight at room temperature</li> </ul>
<b>Extraction</b>	<ul style="list-style-type: none"> <li>Add to separatory funnel along with 30 mL of hexane containing 12.5 mg/mL BHT and shake to extract</li> <li>Wash hexane fraction with 20 mL of water:50% w/v KOH (85:15)</li> </ul>
<b>Sample Preparation</b>	<ul style="list-style-type: none"> <li>Remove 10 mL of hexane extract and dry under nitrogen</li> <li>Reconstitute in 1 mL of acetonitrile:water (70:30) with 5-minute sonication</li> <li>Filter through 0.45 µm PTFE membrane</li> </ul>

**Table 2. Found Concentrations (IU/mL) for Vitamin D with % RSD (n=3)**

Sample	Vitamin D <sub>2</sub>	Vitamin D <sub>3</sub>
Whole milk	0.16 (41%)	4.76 (4%)
Infant formula	0	5.29 (0.2%)

## Results

Chromatography of the vitamin D and hydroxyvitamin D compounds is presented in Figures 2 and 3, respectively. Analyses of vitamin D<sub>2</sub> and D<sub>3</sub> in extracted milk and infant formula are shown in Figures 4 and 5, respectively. No Vitamin D<sub>2</sub> was found in infant formula while a low amount of D<sub>2</sub> was found in milk. Concentrations of Vitamin D<sub>3</sub> in milk and infant formula were found to be close to the expected fortification levels of at least 4.2 IU/mL (Table 2). No hydroxyvitamin D forms were found in samples of either milk or infant formula. The lowest level for calibration curve was 0.2 IU/mL which would correspond to 0.5 ng of vitamins per mL of milk.

## Conclusions

Presented here is a method for the determination of nutritionally important vitamin D forms as standards and extracted from milk and infant formula after saponification following AOAC method 2011.11.

**Figure 2. UHPLC/MS Analysis Vitamin D Standards at 1 IU/mL on Titan C18**

Vitamin D<sub>3</sub> is shown as 385/259 parent/daughter transition and vitamin D<sub>2</sub> is shown as 397/125 transition using the same Y-scale.

column: Titan C18, 10 cm × 2.1 mm I.D., 1.9 μm (577124-U)

mobile phase: [A] 0.1% (v/v) formic acid in methanol: water (20:80);  
[B] 0.1% (v/v) formic acid in methanol

gradient: 60 to 90% B in 0.4 min, to 100% B in 0.3 min, held for 7.8 min,  
to 60% B, held for 1.5 min


flow rate: 0.25 mL/min for 5.55 min, to 0.50 mL/min for 3.84 min,  
to 0.25 mL/min for 0.7 min

column temp.: 25 °C


detector: MS/MS, APCI(+), MRM as indicated in Table 1

injection: 5 μL


sample: vitamin D<sub>2</sub> and D<sub>3</sub>, 1 IU/mL

**Figure 3. UHPLC/MS Analysis of 25-OH-Vitamin D Standards at 250 ng/mL on Titan C18**


25-OH-Vitamin D<sub>3</sub> (red trace) is shown as 383/257 parent/daughter transition and 25-OH-vitamin D<sub>2</sub> (blue trace) is shown as 395/269 transition using the same Y-scale. Conditions same as in Figure 2.

**Figure 4. UHPLC/MS Analysis Vitamin D in Milk on Titan C18**

Vitamin D<sub>3</sub> is shown as 385/259 parent/daughter transition and vitamin D<sub>2</sub> is shown as 397/125 transition using the same Y-scale. Conditions as in Figure 2.

**Figure 5. UHPLC/MS Analysis Vitamin D in Infant Formula on Titan C18**

Vitamin D<sub>3</sub> is shown as 385/259 parent/daughter transition and vitamin D<sub>2</sub> is shown as 397/125 transition using the same Y-scale. Conditions as in Figure 2.



### + Featured Products

Description	Cat. No.
<b>UHPLC Column</b>	
Titan C18, 10 cm × 2.1 mm I.D., 1.9 μm particle size	577124-U
<b>Mobile Phase Components</b>	
Methanol, LC-MS Ultra CHROMASOLV®, tested for UHPLC-MS, 1 L, 2 L	14262
Water, LC-MS Ultra CHROMASOLV, tested for UHPLC-MS, 1 L, 2 L	14263
Formic acid, LC-MS Ultra eluent additive, 1 mL, 2 mL	14265
<b>Standards</b>	
Ergocalciferol (vitamin D <sub>2</sub> )	47768
Cholecalciferol (vitamin D <sub>3</sub> )	47763
Vitamin D <sub>3</sub> (6,19,19-d <sub>3</sub> ) solution	740284
Vitamin D <sub>2</sub> (6,19,19-d <sub>3</sub> ) solution	739847
25-Hydroxyvitamin D <sub>2</sub>	17937
25-Hydroxyvitamin D <sub>3</sub>	17938
<b>Saponification and Extraction Solvents and Reagents</b>	
Hexane, CHROMASOLV, for HPLC, ≥97%	34859
Water, CHROMASOLV Plus, for HPLC, ≥99.9%	34877
Acetonitrile, CHROMASOLV Plus, for HPLC, ≥99.9%	34998
Reagent alcohol (ethanol)	362808
Pyrogallol, ACS reagent	254002
Potassium hydroxide, ACS reagent, ≥85%, pellets	221473
2,6-Di-tert-butyl-4-methylphenol (BHT), ≥99.0% (GC), powder	B1378
<b>Glassware and Accessories</b>	
Aldrich® separatory funnel with glass stopcock, capacity 250 mL	Z548030
Acrodisc® syringe filters; PTFE membrane, diam. 13 mm, pore size 0.45 μm	Z259888