


# Bromate with 887 Professional UV/VIS Detector and 886 Professional Thermostat / Reactor according to ISO 11206



Bromate is a disinfection byproduct generated by ozonation during drinking water purification. It is undesirable, because it is a suspected human carcinogen. ISO 11206 describes its determination by IC with UV detection after post-column reaction. Bromate oxidizes iodide in acid solution to iodine. This forms with iodide the triiodide ( $I_3^-$ ) which absorbs at 352 nm.

## Results

Mineral water spiked

Bromate

Spike: 0.5  $\mu\text{g/L}$

0.501  $\mu\text{g/L}$

# Method description

## Sample

Drinking water, mineral water

## Sample preparation

Direct injection

## Column

Phenomenex Star-Ion™ A300 6.1005.100  
HC - 100/10.0

Metrosep BP 1 Guard/2.0 6.1015.000

## Solutions

Eluent 100 mmol/L H<sub>2</sub>SO<sub>4</sub> /  
19.3 µmol/L ammonium  
heptamolybdate

Post-column reagent 0.27 mol/L KI

## Analysis

UV detection 352 nm

## Parameters

Flow rate column 0.8 mL/min

Flow rate PCR 0.2 mL/min

Injection volume 1000 µL

P<sub>max</sub> 20.0 MPa

Recording time 16 min

Column temperature 30 °C

PCR temperature 30 °C

Measuring duration 300 ms

## Instrumentation

Professional IC Cation – HP Gradient 2.850.1220



858 Professional Sample Processor 2.858.0020

887 Professional UV/VIS Detector 2.887.0010

886 Professional Reactor 2.886.0110

## Calibration curve

Standard Number	Concentration [µg/L]
1	5.0
2	4.5
3	4.0
4	3.5
5	3.0
6	2.5
7	2.0
8	1.5
9	1.0
10	0.5
11	0.1



m