

REACH Polymer Status determined with GPC/SEC

Application Note Chemical Manufacturing

Authors

Daniela Held and Friedhelm Gores contact: <u>dheld@pss-polymer.com</u> and <u>fgores@pss-polymer.com</u>

Introduction

REACH (short for **R**egistration, **E**valuation, **A**uthorisation and restriction of **Ch**emicals) came into force on June 1st 2007 and applies in all EU member states. Polymers are "special substances" in terms of REACH. However, companies have to prove

Polymers are "special substances" in terms of REACH. However, companies have to prove that their products/educts are polymers.

A polymer molecule in the sense of REACH is a molecule that contains a sequence of at least 3 monomer units M, covalently bound to at least one other monomer unit or reactant (initiator). This definition is the so called (3M+1) – rule.

Additional requirements are:

a) Over 50 percent of the weight for that substance consists of polymer molecules (> 3M+1).b) The amount of polymer molecules with the same chain length must be less than 50 weight percent of the substance.

The preferred method to identify whether a substance falls under the definition of a polymer is Gel Permeation Chromatography (GPC), also known as Size Exclusion Chromatography (SEC). GPC/SEC separates substances by chain length enabling percent weight per chain length to be determined. One result of a GPC/SEC experiment is a table with molar mass/number of monomer units and peak area. Table 1, with 2 examples for substances which are not polymers and one substance which is a polymer, illustrates the REACH criteria.

Monomer units	Sample A	Sample B	Sample C	
M =1	0	35%	0%	
М=2	5%	20%	10%	
М=3	55%	10%	20%	
M=4	20%	10%	30%	
M=5	10%	10%	20%	
М=6	5%	10%	10%	
M=7	5%	5%	5%	
M=8	0%	0%	5%	
Polymer definition	failed (b)	failed (only 45% > 3M+1)	passed	

Table 1: REACH criteria discussed with 3 different examples



PSS Polymer Standards Service GmbH In der Dalheimer Wiese 5 55120 Mainz | Germany
 Phone
 +49 6131 96239-0

 Fax
 +49 6131 96239-11

 E-Mail
 info@pss-polymer.com

 Web
 www.pss-polymer.com

Polymer Standards Service-USA, Inc. 160 Old Farm Rd, Suite A Amherst | MA 01002 | USA
 Phone
 +1 413 835-0265

 Fax
 +1 413 835-0354

 E-Mail
 usa@pss-polymer.com

 Web
 www.pss-polymer.com



System Requirements

	Conditions		
Pump	 PSS SECcurity GPC1260 isocratic pump flow rate [mL/min]: 1 mobile phase: THF (other solvents possible) 		
Injection system	PSS SECcurity GPC1260 Autosampler injection volume 100µL 		
Columns	 PSS SDV precolumn (8*50mm) PSS SDV 3µ 1 000 Å (8*300mm) PSS SDV 3µ 1 000Å (8*300mm) PSS SDV 3µ 1 000Å (8*300mm) 		
Calibration	PSS Polystyrene ReadyCal Low		
Loading	5 mg/mL, 20μL injection volume		
Detectors	Refractive index PSS SECcurity 1260 RI		
Software	PSS WinGPC UniChrom		



Procedure, Results & Discussion

Polyols are alcohols with multiple hydroxyl groups. They are used for very different applications, such as sugar substitutes in Feed & Food or as reactants to produce high quality products in polymer chemistry. In order to determine if the above stated REACH polymer requirements for a polyol are true, a GPC/SEC experiment using high resolution SDV columns was performed.

Figure 1 shows the graphical results for a Polyol showing the size-based separation into individual $R-(M)_x$ peaks.



Figure 1: GPC/SEC separation of a polyol into the individual oligomer chains.



PSS Polymer Standards Service GmbH In der Dalheimer Wiese 5 55120 Mainz | Germany
 Phone
 +49 6131 96239-0

 Fax
 +49 6131 96239-11

 E-Mail
 info@pss-polymer.com

 Web
 www.pss-polymer.com

Polymer Standards Service-USA, Inc. 160 Old Farm Rd, Suite A Amherst | MA 01002 | USA
 Phone
 +1 413 835-0265

 Fax
 +1 413 835-0354

 E-Mail
 usa@pss-polymer.com

 Web
 www.pss-polymer.com



The REACH criteria are then applied to all identified individual chains of the Polyol. Table 2 lists peak area (obtained from the detector signal) and the molar mass (obtained using a previously established calibration curve) and comments this with respect to the REACH requirements.

		Content [%]	Relative Molar Mass M _n [Da]*	Comment	Polymer
R-M,	Peak H	0.20	350	Less than 50% below (3M+1)	Yes
R-M-M	Peak G	13.4	650		
R-M-M-M	Peak F	15.7	960		
	Peak E	14.8	1 270		Mara
	Peak D	13.3	1 580	distribution, (3M+1)-Rule	Yes
	Peak C	11.0	1 890		165
	Peak B	9.30	2 200		
	Peak A	22.3	2 550		

Table 2: Content and molar mass for a Polyol and comments to the REACH requirements.

* Polystyrene equivalents

As a

- molar mass distribution with higher molar masses than 3M+1 was measured
- where no individual chain contributes to more than 50%
- and as the combined content of Peak G and H is below 50%

it can be concluded that this Polyol is a polymer according to the definition by REACH.

While it is easy to assign the peaks in the example above, there are also examples where it is more difficult to identify M of the repetition units. In these cases GPC/SEC hyphenated with mass spectrometry, MS, can be a good tool to determine the molar mass of M to ensure a correct REACH data analysis.



PSS Polymer Standards Service GmbH In der Dalheimer Wiese 5 55120 Mainz | Germany
 Phone
 +49 6131 96239-0

 Fax
 +49 6131 96239-11

 E-Mail
 info@pss-polymer.com

 Web
 www.pss-polymer.com

Polymer Standards Service-USA, Inc. 160 Old Farm Rd, Suite A Amherst | MA 01002 | USA
 Phone
 +1
 413
 835-0265

 Fax
 +1
 413
 835-0354

 E-Mail
 usa@pss-polymer.com

 Web
 www.pss-polymer.com