

Application News

No. SCA_190_047

High Performance Liquid Chromatography

Analysis of Capsaicinoids in Chili Products

Introduction

Pungency of chili products depends on the amount of capsaicinoids which are naturally present in bell peppers or chili peppers. The two main components, capsaicin (69%) and dihydrocapsaicin (22%), are almost twice as strong as the capsaicinoids nordihydrocapsaicin (7%), homodihydrocapsaicin (1%) and homocapsaicin (1%), which are smaller in comparison [1]. Therefore, only capsaicin and dihydrocapsaicin were studied to determine capsaicin levels in chili products.




Fig.1 Chemical Structures of Capsaicinoids

The pungency level of a pepper is measured in Scoville Units, defined by Wilbur Scoville. Scoville used a group of tasters to determine the dilution needed to remove the pungency from a sauce or food [2]. Pure capsaicin is rated at 16,000,000 Scoville units.

Method

An isocratic system equipped with a photodiode array detector (PDA) and fluorescence detector was used for the analysis. The method parameters are shown in table 1.

Table 1: Analytical Conditions

System	Nexera X2 + SPD-M30A + RF-20Axs
Column	Shim-Pack GIST C18 2 μ m; 2.1 x 100 mm
Mobil Phase	A: 1 % Acetic Acid in H ₂ O; B: Acetonitrile
Method	A: 60 Vol.-% B: 40 Vol.-%
Flow Rate	0.9 ml/min
Detection	PDA 280 nm RF-20Axs: Ex 280 nm, Em 325 nm
Temp.	50 °C
Inj. Vol.	1 μ l
Run Time	6 min

Fluorescence detection exhibits approximately 16-times higher sensitivity than detection with the PDA. Therefore, the focus of this study was only on the fluorescence detector for further method development.




Fig. 2: Comparison PDA vs. Fluorescence Detector

■ Calibration

The two capsaicinoids capsaicin and dihydrocapsaicin were analyzed. Figure 3 shows a chromatogram of a capsaicinoids standard solution with a concentration of 5 µg/ml of each analyte.




Fig. 3: Chromatogram of a Standard (5 µg/ml)

In order to enable quantification of the capsaicin content, six point calibration curves in a concentration range of 5 – 100 µg/mL were created for each analyte. Calibration curves and the corresponding coefficient of determination R^2 are shown in figure 4.




Figure 4: Calibration Curves of the two Capsaicinoids

■ Sample Preparation and Measurement

Figure 5 shows the sample preparation for the chili products.




Fig. 5: Sample Preparation

An example of a chromatogram of a chili sauce is shown in figure 6. According to literature [3] it is assumed based on the elution order that the peak at 2.6 min is nordihydrocapsaicin.




Fig. 6: Chromatogram of a Chili Sauce

■ Summary

This application note describes the analysis of the two main capsaicinoids of chili peppers using a fluorescence detector in combination with a Nexera X2 UHPLC system. All capsaicinoids elute within 4.5 minutes and were baseline separated. The method was calibrated in a concentration range of 5 to 100 µg/ml and linearity in this range was shown by a coefficient of determination R^2 of ≥ 0.998 for all analytes.

- [1] Bennett DJ, Kirby GW (1968). J. Chem. Soc. C: 442
- [2] Collins MD, Wasmund LM, Bosland PW (1995) HortScience. 30 (1): 137–139
- [3] AOAC Official Method 995.03